

Werkzeuge/ Tools/ Outils

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## Safety instructions

- No liquid may get into the drive unit of the dispenser.
- The work has to be interrupted immediately if you notice unusual noise and/ or emission of liquid from the upper rinsing drilling. The seals must be controlled, and replaced if necessary.
- While working with the dispenser, the user must wear his personal protective equipment, in particular safety goggles.
- Follow the health and safety regulations accident prevention regulations.
- Note the operating instructions of the dispenser unit.
- Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.

### Product information

- The shaft bearing of the dispersion tool consists of a two ball bearings and a single-thrust slide ring seal (SC).
- Content of package: shaft and dispersion head
- The plastic materials used are approved by FDA (Food and Drug Administration Approval).

### Application instructions

- Attention! Never run the dispersion tool dry. Without cooling by the medium the bearing will be destroyed.
- The dispersion tool is not suitable for continuous operation.
- Dispensing causes the medium to heat up.
- Likewise the optimal dispersion duration and rotating frequency must be determined by attempts.
- The function of the dispersion tools depends on the condition of the sharp edges of rotor and stator. Abrasive media can round off these edges fast, whereby the disperse effect diminishes.

Designation	Application	Max. permissible number of revolutions [1/min]	Working range pressure [bar]	Working range vacuum [mbar]	Ultimate fineness emulsions
S 25 KV-18 G	Wet crushing, suspension, emulsifying, disagglomeration (batch-operation)	24000	6	1	1-10 µm
S 25 KV 25 G	Wet crushing, suspension, emulsifying, disagglomeration (batch-operation)	24000	6	1	1-10 µm
S 25 KV25 F	Wet crushing, suspension, emulsifying, disagglomeration (batch-operation)	24000	6	1	1-5 µm

### Warranty

In accordance with **IKA** warranty conditions, the warranty periods 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

### Spare parts list

<b>Item</b>	<b>Designation</b>	<b>Item</b>	<b>Designation</b>
1	Cylindrical pin	12	O-Ring 9,25 x 1,78
4	Deep-groove ball bearing	3003	Pressure spring, right
6	Seeger circlip ring	3004	Thrust collar
8	Stator	3005	Slide ring
9	Rotor	2001	Rotor shaft, compl.
10	O-Ring 4,7 x 1,42	2002	Friction ring seal, compl.
		2003	Shaft tube, compl.

Spare parts diagram see front page.

For orders of spare parts please specify the serial number on the shaft tube, the type of dispersion tool and the designation of the spare part.

### Cleaning, disinfection and sterilization

Allowed procedures	Methods of sterilisation	Notice
Chemical procedures	By germ - killing disinfectants alcohol, phenol,formalin...	For ball bearings <b>not</b> allowed. Note the chemical resistance of ball bearings and o-rings.
Hot air	Disinfectant remainders are to be removed with germ-free water.	Life span of the bearings shortens

For pre-cleaning the dispersion tool, operate it in a dissolving for dissolving the residues of the substance. Due to the large flow rate, rotor and stator are cleaned to a large extent. The dispersion tool must be divided and cleaned immediately after working, so that the adhering residues of the substance don't cultivate unwanted bacterial cultures.

For cleaning purposes disassemble the dispersion tool as described below (items see front page):

- Unscrew the rotor (9) from the shaft (2001) using the generator and shaft key included in the tool set. The rotor is held fast by means of the generator key and screwed off from above by turning the shaft using the shaft key.
- Unscrew the stator (8) from the shaft (2) (left-hand thread). Insert transverse rib of the shaft key into the front-end slot of the shaft tube. By using the generator key the stator may now be screwed off in a clockwise direction. In doing so, the dispersion tool and the stator should be held up in vertical position, thus preventing the spring of the slide ring sealing from falling out.
- Tap the shaft (2001) and the ball bearing (4) downward in the direction of the stator thread out of the shaft tube.
- In this position the parts of the slide ring seal can be removed from the shaft (2001).
- The slide rings and O-rings must be replaced in the case of wear. The sliding surfaces of the slide and counter ring must be clean.

The assembly of the dispersion tool takes place in reverse order. Make certain that the O-ring (10) and the thrust collar (3004) are properly fitted during assembly.

Attention! The rotor must not be tightened too strongly, because otherwise the thread will be damaged.

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