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Dissolved oxygen probe StirrOx G with integrated stirrer

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
**Construction of the dissolved
oxygen probe StirrOx G**

- | | | |
|---|---|--|
| ① | Membrankopf WP-ST | <i>Membrane head WP-ST</i> |
| ② | Temperaturfühler | <i>Temperature probe</i> |
| ③ | Schaftgehäuse | <i>Shaft housing</i> |
| ④ | Abschlusskopf | <i>Closing head</i> |
| ⑤ | Arbeitselektrode (Goldkathode) | <i>Working electrode (gold cathode)</i> |
| ⑥ | Gegenelektrode (Bleianode) | <i>Counter electrode (lead anode)</i> |
| ⑦ | Isolator | <i>Isolator</i> |
| ⑧ | Ein/Aus-Taster | <i>On/off key</i> |
| ⑨ | Rührpaddel RP-ST | <i>Stirring paddle RP-ST</i> |
| ⑩ | Stecker 8-polig | <i>8 pin plug</i> |
| ⑪ | Stecker 2-polig
(Stromversorgung für Rührer) | <i>2 pin plug
(power supply for stirrer)</i> |

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The membrane cap that is mounted on the sensor for delivery serves mainly as a transport protection. Depending on the duration of the transport and storage period, it may have a shortened operational lifetime.

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*Connect probe to the meter.
The probe is immediately ready for measurement.*

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A polarization of the probe is not required.

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For calibration please refer to the operation manual of the meter.

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*Note:
Do not switch on the stirrer when the probe is in the calibration beaker. The stirrer cannot rotate there.*

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For the necessary minimum flow, the probe has been provided with a built-in stirring system.

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*To turn on the stirrer, press the on/off key.
To turn it off, press the on/off key again.*

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Recommendation:

After each series of measurements (before storing the probe) clean the probe (exterior cleaning).

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Important note:

The StirrOx G has a displacement volume of 5.6 ml. Therefore, the sample bottle is closed free of air bubbles after inserting the stopper.

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If other bottles than the Karlsruhe bottles available at WTW (see chapter "Recommended accessories") are used, please use only bottles with standard grinding NS 19. The stopper must have a displacement volume between 5.6 and 8.6 ml.

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Cleaning procedure

Lime

Immerse for 1 minute into acetic acid w = 25 %.

Grease/Oil

Rinse with warm water and household cleaning solution.

Rinse thoroughly with deionized water.

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Store the probe in the calibration beaker.

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Humidity: Keep air in calibration beaker moist

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See application report 895234 (in the appendix of this manual):
"Regeneration of galvanic oxygen sensors".

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Error symptom	Possible cause	Elimination
<i>Instrument displays 0.0 mg/l or 0 % Sat O₂ (probe in air)</i>	<ul style="list-style-type: none"> ● <i>No connection between meter and probe</i> ● <i>Defective cable</i> 	<p><i>Control connection between meter and probe</i></p> <p><i>Send probe to WTW</i></p>
<i>Probe cannot be calibrated</i>	<ul style="list-style-type: none"> ● <i>Contaminated membrane head</i> 	<i>Refill probe according to operation manual, wait 60 min. and recalibrate</i>
<i>Probe cannot be calibrated even after exchange of electrolyte and membrane head</i>	<ul style="list-style-type: none"> ● <i>Electrode contaminated or probe poisoned</i> 	<i>Cleaning of electrode</i>
<i>Leakage is displayed</i>	<ul style="list-style-type: none"> ● <i>Membrane head has not been really tightened</i> ● <i>Hole in the membrane</i> 	<p><i>Screw membrane head tight</i></p> <p><i>Exchange membrane head according to operation manual</i></p>
<i>Wrong temperature indication</i>	<ul style="list-style-type: none"> ● <i>Defective temperature probe</i> 	<i>Send probe to WTW</i>
<i>Mechanical damage of the probe</i>		<i>Send probe to WTW</i>
<i>Instrument displays OFL</i>	<ul style="list-style-type: none"> ● <i>Display range exceeded</i> ● <i>Short circuit between gold cathode and lead anode</i> 	<i>Regenerate probe; if meter still displays OFL send probe to WTW</i>
<i>Stirrer does not rotate</i>	<ul style="list-style-type: none"> ● <i>Power supply has not been turned on</i> ● <i>Stirring system defective</i> 	<p><i>Turn on power supply</i></p> <p><i>Send probe to WTW</i></p>

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I_0 Static current

I_1 Current during normal operation

I_2 Current when the stirrer is blocked

U_e / V	I_0 / mA	I_1 / mA	I_2 / mA
5,0	< 4,0	< 38	< 189
6,0	< 4,8	< 35	< 147
7,0	< 5,8	< 32	< 137
8,0	< 6,7	< 29	< 121
9,0	< 7,7	< 28	< 108
10,0	< 8,6	< 27	< 100
11,0	< 9,5	< 27	< 95
12,0	< 10,4	< 27	< 89

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Regeneration of galvanic oxygen sensors

WTW delivers the sensor ready to use.

An exchange of the membrane and the electrolyte solution is necessary in case of:

- damaged membrane,
- strongly contaminated membrane
- used up electrolyte solution.

Instrum. and acc.

Exchange membrane heads
Cathode cleaner SF 300
Paper towels
Beakers

Reagents

Deionized water
Cleaning solution RL/G
Electrolyte solution ELY/G

Hints

For all maintenance operation disconnect the sensor from the instrument.

The sensor is ready for operation after regeneration in approx. 30 to 50 minutes.

For measurements in agents with very low oxygen concentrations (< 0,5 % sat), rest the sensor approx. 12 hours.

The electrolyte solution ELY/G is strongly alkaline!


- Keep out of reach of children!
- Wear suitable gloves and eye/face protection!
- After contact with eyes and skin wash immediately with plenty of water!
- In case of contact with eyes seek medical advice!

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
Tel.: 0881/183-0 • Fax (0881) 6 25 39


Regeneration of the oxygen sensor


Disconnect plugs from the instrument.


Unscrew membrane head.


Rinse sensor head with deionized water.


Polish impurities from the cathode with cathode cleaner (wet grinding foil) under low pressure. Don't use grinding paper or glass fibre pen!

Rinse sensor head with deionized water.


Wipe off counter electrode with fuzzy-free paper towel and remove carefully white coating.


Immerse sensor head including the counter electrode into cleaning solution RL/G for 1 to 3 minutes.


Please note the safety guidelines on the bottle!


Rinse sensor head several times with deionized water.


Shake carefully off water drops.


Fill a new membrane head carefully with electrolyte solution.


Remove existing air bubbles by careful knocking.

Rinse sensor head with electrolyte solution.

- Screw membrane head on the shaft.
- Hold sensor inclined. Excessed electrolyte solution is pushed out.
- Screw on membrane head fingertight with a paper towel.

Small air bubbles do not interfere.

Refilling only in case of:

- Big air bubbles,
- Air bubbles at the gold cathode.

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Wissenschaftlich-Technische Werkstätten GmbH & Co. KG