



## E 136







Insert  
for optimum loading of 56 half petri dishes.


- Suitable for 56 half petri dishes with Ø 100 mm
- equipped with 56 holders, height 70 mm, distance apart approx. 26 mm
- For use in lower basket



## Product details - E 136

EAN: 4002511115941 / Article no. 69513601D / Mat. no. 03830280

Product category 	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine/accessory categorisation 	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category 	
Insert/module for petri dishes	•
Product properties 	
Material	Stainless steel
Colour	Stainless steel
Distance between rails in mm	26
Rail height in mm	70
Area of application 	
Reprocessing laboratory glassware	•
Capacity 	
Petri dishes [number]	56

Dimensions and weight 

External dimensions, height in mm	485
External dimensions, width in mm	145
External dimensions, depth in mm	445
Net weight in kg	2,100

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432

Fax: (01) 4523967

E-mail: [info@labunlimited.com](mailto:info@labunlimited.com)

Web: [www.labunlimited.com](http://www.labunlimited.com)

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30

Fax: 08452 30 50 30

E-mail: [info@labunlimited.co.uk](mailto:info@labunlimited.co.uk)

Web: [www.labunlimited.co.uk](http://www.labunlimited.co.uk)