

# PerkinElmer/Sciex

Elan 250/500/5000/6000/6100/9000/DRC

✓ Standard Configuration • Optional Configuration



## CONES

Part #	Description	Cross Ref.	5000	6X00/9000	DRC
SC2003-Ni	Ni Sampler Cone	WE013087	✓		
SC2004-Ni	Ni Skimmer Cone	N8122000	✓		
SC2007-Pt	Pt Sampler Cone	N8122025	✓		
SC2010-Pt	Pt Skimmer Cone	N8122035	✓		
SC2011-Al	Al Sampler Cone	-		✓	✓
SC2011-Ni	Ni Sampler Cone	WE021140		✓	✓
SC2012-Al	Al Skimmer Cone	-		✓	✓
SC2012-Ni	Ni Skimmer Cone	WE021137		✓	✓
SC2013-Pt	Pt Sampler Cone	WE027802		✓	✓
SC2013-Pte	e • SPEC: Pt Sampler Cone	WE027802		✓	✓
SC2014-Pt	Pt Skimmer Cone	WE027803		✓	✓
SC2014-Pte	e • SPEC: Pt Skimmer Cone	WE027803		✓	✓
SC3000-Pte	Pt Cone Kit (SC2013-Pte, SC2014-Pte, O-Rings)*	-		✓	✓
SC5011	Sampler Cone O-rings (pkg. 5)	N8120511		✓	✓
SC5012	Skimmer Cone O-rings (pkg. 5)	N8120512		✓	✓
SC6000	Cone removal / insertion tool	WE017142		✓	✓

\*Special pricing - Available only on the web - [www.SpectronUS.com](http://www.SpectronUS.com)

## NEBULIZERS

Part #	Description	Cross Ref.	5000	6X00/9000	DRC
AR30-07-FC3E	Conikal Type K (3mL/min), Glass	N0681574	✓	✓	
500-80Q	Concentric Type K (3mL/min), Quartz	N0681574	✓	✓	
500-72Q	Concentric Type C (3mL/min), Quartz	WE024371	✓	✓	✓
AR30-1-FC3E	Conikal Type C (3mL/min), glass	N8102011	✓	✓	✓
ES2040-xx	Concentric PFA-ST Neb, need ES2542	N8122192	✓	✓	✓
ES2020	Concentric PFA-20 Nebulizer	N8122350	✓	✓	✓
ES2002	Concentric PFA-100 Nebulizer	N8122351	✓	✓	✓
ES2542	Teflon union for PFA-ST, & autosampler	N8122355			

# PerkinElmer/Sciex

Elan 250/500/5000/6000/6100/9000/DRC

✓ Standard Configuration • Optional Configuration

## RF COILS

Part #	Description	Cross Ref.	5000	6X00/9000	DRC
70-900-2005C	RF Coil Assembly, Copper	WE012411	✓		
70-900-2006C	RF Coil Assembly, Copper	WE021816		✓	✓



## MULTIPLIERS

Part #	Description	Cross Ref.	5000	6X00/9000	DRC
AF210	ETP Electron Multiplier, Dual Mode	N8125001		✓	✓
AF217 <sup>+</sup>	ETP Electron Multiplier, Dual (20% longer lifetime)	N8125050		✓	✓
AF570	ETP Electron Multiplier, 5000	N8101069	✓		
AF571	ETP Electron Multiplier, 5000a	N8122189	✓		



\* Instruments with serial numbers that do not begin with an "A" (before 2005) require upgrade W101-3195 before installing this detector

## TORCHES, SPRAY CHAMBERS & ACCESSORIES

Part #	Description	Cross Ref.	5000	6X00/9000	DRC
400-10	Standard Quartz Torch	N8122006	✓	✓	✓
400-20J	Quartz Jacketed Cyclonic Chamber (water-cooled)	N8120181			✓
400-20	Quartz Cyclonic Chamber, w/ ball joint	WE025221			✓
400-36	Torch O-ring Kit (pkg 5)	N8120100	✓	✓	✓
400-37	Injector Support, Elan DRC	WE023951			✓
400-35	Injector Support, Torch Base	N8120116	✓	✓	
400-30	Qtz Injector, ball joint 2.0mm (standard)	WE023948			✓
400-32	Qtz Injector, ball joint 1.5mm (low flow)	WE027005			✓
400-31	Qtz Injector, ball joint 0.85mm (volatile organic solv)	WE027030			✓
320-35	Qtz Injector, 2.0mm i.d. (clinical samples)	N8125029	✓	✓	
320-34	Qtz Injector, 1.6mm i.d. (organic solvents)	00473292	✓	✓	
320-33	Qtz Injector, 1.2mm i.d. (organic solvents)	N0681631	✓	✓	
320-53	Alumina Injector, 2.0mm (aqueous solutions)	N8126041	✓	✓	
320-52	Alumina Injector, 1.5mm (low flow)	N8126040	✓	✓	
320-51	Alumina Injector, 0.85mm (volatile organic solv)	N8126039	✓	✓	
20-809-0297	Glass Jacketed Cyclonic Chamber (water-cooled)	N8120181			✓
20-809-0224	Glass Cyclonic Chamber, 50mL	N8120150	✓	✓	
ES-1124-0201	Sapphire Injector, 2.0mm (specify instrument)	N0695495	✓	✓	✓



Distributed by:



**Carl Stuart Limited**

ADVANCED APPLIED TECHNOLOGIES

Contact Us:

Int Ph: 01 4523432

UK Ph: 08452 30 40 30

Web: [www.carlstuart.com](http://www.carlstuart.com)

Email: [info@carlstuart.com](mailto:info@carlstuart.com)