


E 110

Insert
for the optimum loading of beakers from 250 to
600 ml.



- Suitable for 10 glass beakers between 250 - 600 ml
- equipped with 10 x 3 stay bars
- For use in lower basket



Product details - E 110

EAN: 4002511115736 / Article no. 69511001D / Mat. no. 03808390

Product category		^
Washer-disinfectors, laboratory		•
Large-capacity washer-disinfectors, laboratory		•
Machine/accessory categorisation		^
G 7825		•
PG 8527		•
PG 8535		•
PG 8536		•
PG 8583		•
PG 8583 CD		•
PG 8593		•
Product category		^
Insert/module for beakers		•
Product properties		^
Material	Stainless steel	
Colour	Stainless steel	
Area of application		^
Reprocessing laboratory glassware		•
Dimensions and weight		^
External dimensions, height in mm	460	
External dimensions, width in mm	230	
External dimensions, depth in mm	175	

Net weight in kg

0,970

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432

Fax: (01) 4523967

E-mail: info@labunlimited.com

Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30

Fax: 08452 30 50 30

E-mail: info@labunlimited.co.uk

Web: www.labunlimited.co.uk