


## E 144

## Insert







For the optimum loading of up to 18 beakers.

- Suitable for 18 glass beakers
- equipped with 18 x 3 stay bars
- For use in upper and lower basket



## Product details - E 144

EAN: 4002511115996 / Article no. 69514401D / Mat. no. 03808710

Product category 	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine/accessory categorisation 	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category 	
Insert/module for beakers	•
Product properties 	
Material	Stainless steel
Colour	Stainless steel
Area of application 	
Reprocessing laboratory glassware	•
Dimensions and weight 	
External dimensions, height in mm	131
External dimensions, width in mm	200
External dimensions, depth in mm	445

Net weight in kg

1,370

Subject to technical changes; no liability accepted for the accuracy of the information given!

Distributed by:

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432  
Fax: (01) 4523967  
E-mail: [info@labunlimited.com](mailto:info@labunlimited.com)  
Web: [www.labunlimited.com](http://www.labunlimited.com)

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30  
Fax: 08452 30 50 30  
E-mail: [info@labunlimited.co.uk](mailto:info@labunlimited.co.uk)  
Web: [www.labunlimited.co.uk](http://www.labunlimited.co.uk)